

SUBJECT	CONVERSION	SHEET #	CONVERSION.FIE
---------	------------	---------	----------------

Conversion chart showing approximate equivalents of metric and gauge/inch sizes

Ø		Length	
mm	Gauge	mm	Inches
0.30	30G	3	1/8"
0.35	28G	4	1/6"
0.40	27G	5	1/5"
0.45	26G	6	1/4"
0.50	25G	8	5/16"
0.55	24G	10	3/8"
0.60	23G	13	1/2"
0.70	22G	15	5/8"
0.80	21G	20	3/4"
0.90	20G	22	7/8"
1.00 / 1.10	19G	25	1"
1.20 / 1.30	18G	30	1 1/4"
1.40 / 1.50	17G	35	1 3/8"
1.60	16G	40	1 1/2"
1.80	15G	45	1 3/4"
2.00 / 2.10	14G	50	2"
2.40	13G	60	2 1/2"
2.80	12G	70	2 3/4"
3.00	11G	75	3"
3.50	10G	80	3 1/4"
4.00	8G	90	3 1/2"
4.50	7G	100	4"
5.00	6G	110	4 1/4"
		120	4 3/4"
		125	5"
		150	6"
		200	8"

All this information are for reference only. They have no legal or contractual commitment Unimed SA.