







LONGUEURS DES POINTES DE NOS AIGUILLES A INJECTION

POINT LENGTHS OF OUR INJECTION NEEDLES

SPITZEN-LÄNGEN UNSERER INJEKTIONS-KANÜLEN



Ø ext. [mm]	Biseau moyen Regular bevel Normale Spitze		Biseau court Short bevel Kurze Spitze		Biseau extra-court Very short bevel Extrakurze Spitze		Trocart & crayon Trocar & pencil point Trokar & konische	
	B [mm]	α [~]	B [mm]	α [~]	B [mm]	α [~]	C [mm]	β [~]
0.20	0.8	14°	0.6	18°	0.3	34°	0.6	9°
0.25	1.0	14°	0.7	20°	0.4	32°	0.7	10°
0.30	1.2	14°	0.9	18°	0.5	31°	0.9	9°
0.35	1.4	14°	1.0	19°	0.6	30°	1.0	10°
0.40	1.6	14°	1.2	18°	0.7	30°	1.2	9°
0.45	1.8	14°	1.3	19°	0.8	29°	1.3	10°
0.50	2.0	14°	1.5	18°	0.9	29°	1.5	9°
0.55	2.2	14°	1.6	19°	1.0	29°	1.6	10°
0.60	2.4	14°	1.8	18°	1.0	31°	1.8	9°
0.65	2.6	14°	2.0	18°	1.1	31°	2.0	9°
0.70	2.8	14°	2.2	18°	1.2	30°	2.2	9°
0.80	3.2	14°	2.4	18°	1.3	32°	2.4	9°
0.90	3.6	14°	2.7	18°	1.5	31°	2.7	9°
1.00	4.0	14°	3.0	18°	1.7	30°	3.0	9°
1.10	4.4	14°	3.3	18°	1.9	30°	3.3	9°
1.20	4.8	14°	3.3	20°	2.0	31°	3.3	10°
1.30	5.2	14°	3.5	20°	2.2	31°	3.5	11°
1.40	5.6	14°	3.8	20°	2.4	30°	3.8	10°
1.50	5.6	15°	3.8	22°	2.6	30°	3.8	11°



Ø ext. [mm]	Biseau moyen Regular bevel Normale Spitze		Biseau court Short bevel Kurze Spitze		Biseau extra-court Very short bevel Extrakurze Spitze		Trocart & crayon Trocar & pencil point Trokar & konische	
	B [mm]	α [~]	B [mm]	α [~]	B [mm]	α [~]	C [mm]	β [~]
1.60	5.9	15°	4.0	22°	2.6	32°	4.0	11°
1.70	5.9	16°	4.0	23°	2.7	32°	4.0	12°
1.80	6.2	16°	4.2	23°	2.8	33°	4.2	12°
1.90	6.2	17°	4.2	24°	3.0	32°	4.2	13°
2.00	6.5	17°	4.4	24°	3.0	34°	4.4	13°
2.10	6.5	18°	4.7	24°	3.1	34°	4.7	13°
2.20	6.8	18°	4.9	24°	3.2	35°	4.9	13°
2.30	7.1	18°	5.1	24°	3.4	34°	5.1	13°
2.40	7.3	18°	5.3	24°	3.5	34°	5.3	13°
2.50	7.6	18°	5.6	24°	3.7	34°	5.6	13°
2.60	8.0	18°	5.8	24°	3.8	34°	5.8	13°
2.80	8.6	18°	6.2	24°	4.1	34°	6.2	13°
3.00	9.2	18°	6.7	24°	4.4	34°	6.7	13°
3.50	10.0	19°	7.9	24°	5.2	34°	7.9	12°
4.00	11.0	20°	9.0	24°	5.9	34°	9.0	13°
4.50	12.1	20°	9.8	25°	6.5	35°	9.8	13°
5.00	13.2	21°	10.5	25°	7.1	35°	10.5	13°
5.50	14.3	21°	11.2	26°	7.7	36°	11.2	14°
6.00	15.3	21°	11.8	27°	8.3	36°	11.8	14°
6.50	16.2	22°	12.4	28°	8.8	36°	12.4	15°
7.00	17.1	22°	12.9	28°	9.3	37°	12.9	15°
7.50	18.0	23°	13.4	29°	9.8	37°	13.4	16°
8.00	18.8	23°	13.9	30°	10.2	38°	13.9	16°