



CONVERSION CHART

Showing approximate equivalents of metric and gauge/inch sizes

<u>Ø mm</u>	<u>Gauge</u>	<u>Length mm</u>	<u>Inches</u>
0.30	30G	3	1/8"
0.35	28G	6	1/4"
0.40	27G	8	5/16"
0.45	26G	10	3/8"
0.50	25G	13	1/2"
0.55	24G	15	5/8"
0.60	23G	20	3/4"
0.70	22G	22	7/8"
0.80	21G	25	1"
0.90	20G	30	1 1/4"
1.00/1.10	19G	35	1 3/8"
1.20/1.30	18G	40	1 1/2"
1.40/1.50	17G	45	1 3/4"
1.60	16G	50	2"
1.80	15G	60	2 1/2"
2.00/2.10	14G	70	2 3/4"
2.40	13G	75	3"
2.80	12G	80	3 1/4"
3.00	11G	90	3 1/2"
3.50	10G	100	4"
		110	4 1/4"
4.00	8G	120	4 3/4"
4.50	7G	125	5"
5.00	6G	150	6"
		200	8"

Please refer to page 4 of our medical needle catalogue no. 2008A for a list of stock sizes and other standard sizes. Within a diameter-range of 0.30 mm to 5.00 mm practically any length up to 1000 mm can be manufactured on special order (with or without a hub).